

PATGuard 3 Tester Compatibility

The program can **receive** downloaded data from the following testers:

- Seaward Apollo 600
- Seaward Apollo 500
- Seaward PrimeTest 250+
- Seaward PrimeTest 350
- Seaward PrimeTest 300
- Seaward Europa Series
- Seaward Supernova Series
- Seaward IT1000
- Seaward PAT1000 Series
- Seaward PAT 2000 Series
- Seaward SEN60204
- Clare HAL 101/102/103
- Clare SafeCheck 8
- Clare SafeTest Luminaire
- Clare SafeTest Manufacturing
- Ethos 9300 / 9500 / 9600
- MPAT60
- Emona Maestro
- Metrohm MicroPAT+
- Megger PAT4*
- Robin SmartPAT 5000/5500
- Metrel OmegaPAT
- Metrel BETA PAT

****Requires download through a parallel port via a parallel to serial adapter – not available from Seaward – contact Megger for download cable and adaptor requirements.***

The program can **upload** to the following PAT testers:

- Seaward Apollo 600
- Seaward Apollo 500
- Seaward PrimeTest 350
- Seaward PrimeTest 300
- Seaward Europa Series
- Seaward Supernova Series
- Seaward IT1000
- Seaward PAT2000i/iee
- Ethos 9500 / 9600
- Emona Maestro